- Resolution to 5³/₄ digits (1 part out of ±300,000)
- Accuracy to ±0.10% of reading
- Two line display
- Peak capture
- Analog voltage outputs
- IEEE-488 and serial interface

Product Description

Designed for demanding laboratory and materials analysis applications, the Model 450 High Precision Hall Effect Gaussmeter has high resolution and low noise floor. High performance is combined with ease of use – frequently used features can be accessed by one or two keystrokes, the message-based display is bright and the keypad well-defined. As an added advantage, the Model 450 includes a Lake Shore Hall probe.

Measurement Modes

The Model 450 operates in DC, RMS, and Peak modes, with superior accuracy and resolution in DC measurement mode. Measurements to 5³/₄ digits are possible due to the low noise floor. With low noise and high stability, the Model 450 is ideal for field control or mapping applications. Changing fields that are often used in material analysis systems can also be measured up to 18 times per second over computer interface with excellent resolution.

Best suited for linear power supply measurements or measurements of magnets and solenoids driven at line frequency, RMS mode measures periodic AC fields from 10 Hz to 400 Hz. Instrument circuitry accommodates wave forms with crest factors up to 7, with true RMS conversion.

Peak circuitry in the Model 450 captures single event peaks or monitors the peak amplitude of periodic wave forms from 10 Hz to 400 Hz, with reproducible single peak measurements to 5 ms rise time. Instrument software accommodates indefinite hold time with no decay. The Model 480 Fluxmeter is a good choice if faster peak or RMS measurements are required.

Model 450 Gaussmeter



Range and Resolution

When used with appropriate probes, the Model 450 offers full scale ranges from 300 mG to 300 kG. With 5³/₄ digit resolution, DC field variations approaching 0.010 mG can be detected; in larger DC fields, resolution to 1 part in 300,000 is possible. For RMS and peak measurements, resolution is 4³/₄ digits or 1 part in 30,000. The filter feature of the Model 450 improves resolution in noisy environments by taking a running average of field readings. DC mode requires filtering to achieve 5³/₄ digit resolution.

Interface

The Model 450 is equipped with parallel (IEEE-488) and serial (RS-232C) computer interfaces for command and data exchange; maximum reading rate can be achieved with the IEEE-488 interface. Nearly every function of the Model 450 front panel can be performed *via* computer interface. The Model 450 also provides one corrected and one monitor analog voltage output. Corrected for sensor linearity, offset, and temperature effects, the corrected output is a DC voltage proportional to the display reading. The monitor output is a real-time analog voltage proportional to the field.

Probes and Sensors

The Model 450 is compatible with most Lake Shore gaussmeter Hall probes. When ordering the Model 450 with one of the gaussmeter Hall probes on the following page, a discounted package price is available. Lake Shore probes are factory calibrated for accuracy and interchangeability. Factorycalibrated probes feature a programmable read-only memory (PROM) in the probe connector so that calibration data can be read automatically by the instrument. If the probe is equipped with a temperature sensor, the Model 450 reads both temperature and field signal and continuously adjusts the calculated field value. Lake Shore can also custom design a probe to meet your specific application requirements.

Display

The Model 450 has a two line by 20 character vacuum fluorescent display. During normal operation, the display is used to report field readings and give results of other features such as max/min or relative. When setting instrument parameters, the display gives the operator meaningful prompts and feedback to simplify operation. The operator can also control display brightness.

Following are three examples of the various display configurations:

+101.840 kG DC

Normal Reading – the default mode with the display of the live DC field reading.



Max DC Hold On – the maximum value is shown in the lower display while the upper display contains the live DC field reading.



Alarm On – the alarm gives audible and visual indication of when the field value is selectively outside or inside a user specified range. An output relay facilitates pass/fail actuation.

Lake Shore Cryotronics, Inc.

(614) 891-2244

fax: (614) 818-1600


Model 450 Gaussmeter

Model 450 Rear Panel

- Line Input Assembly
- 2 IEEE-488 Interface
- **3** Serial I/O Interface
- 4 Corrected Analog Output
- **5** Monitor Analog Output
- 6 Probe Input



Gaussmeter Hall Probes The Model 450 includes one of the Lake Shore probes listed below – specify probe model number when ordering. See page 24 for details on properly selecting a probe and for a complete listing of available probe models.



¹Temperature compensated probes

30 kG

30 kG

Model 450 Gaussmeter

Model 450 Specifications

General Measurement

Number of inputs: 1

Update rate: 5 rdg/s on display; up to 18 rdg/s with IEEE-488 interface Measurement modes: DC, RMS, Peak Probe compatibility: Standard and custom probes Probe features: Linearity Correction, Temperature Correction, Auto Probe Zero Measurement features: Auto Range, Max Hold, Relative Mode, Filter

Probe connector: 15-pin D style

DC Measurement DC display resolution: 5³/₄ digits with filter, 4³/₄ digits without filter

Probe type Range	Resolution w/ Filter	Resolution w/out Filter
HST Probe		
300 kG	0.001 kG	0.01 kG
30 kG	0.0001 kG	0.001 kG
3 kG	0.00001 kG	0.0001 kG
300 G	0.001 G	0.01 G
HSE Probe		
30 kG	0.0001 kG	0.001 kG
3 kG	0.00001 kG	0.0001 kG
300 G	0.001 G	0.01 G
30 G	0.0001 G	0.001 G
UHS Probe		
30 G	0.0001 G	0.001 G
3 G	0.00001 G	0.0001 G
300 mG	0.001 mG	0.01 mG

DC accuracy: $\pm 0.10\%$ of reading $\pm 0.005\%$ of range

DC temperature coefficient: $\pm 0.05\%$ of reading $\pm 0.003\%$ of range per °C **DC precision:** No filter – ranges $\leq 30 \text{ G} \pm 0.013\%$, ranges $\geq 300 \text{ G} \pm 0.003\%$; Filter¹ – ranges $\leq 30 \text{ G} \pm 0.0025\%$, ranges $\geq 300 \text{ G} \pm 0.0007\%$

AC RMS & Peak Measurement

AC display resolution: 43/4 digits

Probe type Range	RMS Resolution	Peak Resolution
HST Probe		
300 kG	0.01 kG	0.01 kG
30 kG	0.001 kG	0.001 kG
3 kG	0.0001 kG	0.0001 kG
300 G	0.01 G	×
HSE Probe		
30 kG	0.001 kG	0.001 kG
3 kG	0.0001 kG	0.0001 kG
300 G	0.01 G	0.01 G
30 G	0.001 G	×
UHS Probe		
30 G	0.001 G	0.001 G
3 G	0.0001 G	0.0001 G
300 mG	0.01 mG	×

AC RMS frequency range: 10 Hz to 400 Hz

AC RMS accuracy: ±2% of reading (50 Hz to 60 Hz)

AC RMS frequency response: 0 to -3.5% of reading (10 Hz to 400 Hz)

(All AC specifications for sinusoidal input >1% of range)

AC peak accuracy: ±5% typical

AC peak speed: 5 ms

Front Panel

Display type: 2 line \times 20 character, vacuum fluorescent Display resolution: To \pm 5¾ digits Display update rate: 5 rdg/s Displays units: Gauss (G), Tesla (T) Units multipliers: μ , m, k Annunciators: RMS: AC input signal, DC: DC input signal, MAX: Max Hold value, s: Relative reading, R: Remote operation, J: Alarm on Keypad: 21 full travel keys Front panel features: Intuitive operation, display prompts, front panel lockout, brightness control

¹Filter set averaging 8 readings over 1.6 s

Interfaces **RS-232C** capabilities Baud: 300, 1200, 9600 Connector: RJ-11 Update rate: Up to 15 rdg/s **IEEE-488** capabilities Complies with IEEE-488.2 SH1, AH1, SR1, RL1, PP0, DC1, DT0, C0, E1 Software support: LabView™ driver Update rate: Up to 18 rdg/s Alarm Settings: High and low set point, Inside/Outside, Audible Actuators: Display annunciator, beeper Monitor analog output Configuration: Real time analog voltage output Scale: $\pm 3 V = \pm FS$ on selected range Frequency response: DC to 400 Hz Accuracy: Probe dependent **Minimum load resistance:** $1 \text{ k}\Omega$ (short circuit protected) Connector: BNC **Corrected analog output** Configuration: Voltage output generated by DAC Range: ±3 V; ±10 V for the Model 450-10 Scale: User defined Resolution: 0.366 mV Update rate: 5 rdg/s on display; up to 18 rdg/s with IEEE-488 interface Accuracy: ±0.1% of full scale in addition to measurement error Minimum load resistance: 1 kQ (short circuit protected) Connector: BNC

General

Ambient temperature: 15 to 35 °C at rated accuracy; 5 to 40 °C with reduced accuracy Power requirement: 100, 120, 220, 240 VAC (+5%, -10%), 50 or 60 Hz, 20 VA Size: 216 mm W × 89 mm H × 318 mm D (8.5 in × 3.5 in × 12.5 in), half rack Weight: 3 kg (6.6 lb) Approval: CE mark

Ordering Part number 450 450-10	Information Description Model 450 Gaussmeter plus one probe Model 450 Gaussmeter (with corrected analog output set to ± 10 V instead of ± 3 V)
Select a power	configuration
VAC-100	100 VAC, includes U.S. power cord
VAC-120	120 VAC, includes U.S. power cord
VAC-220	220 VAC, includes universal Europe power cord
VAC-240	240 VAC, includes universal Europe power cord
VAC-120-ALL	120 VAC, includes U.S. & universal Europe power cords & all fuses
Accessories in	cluded
115-006	Detachable line cord (U.S.)
115-007	Detachable line cord (European)
4060	Zero gauss chamber
MAN-450	Model 450 user manual
	RJ-11 4-wire cable assembly used with RS-232C interface – cable is 4.3 m (14 ft) long RJ-11 to DB-25 adapter – connects computer to RS-232C port RJ-11 to DE-9 adapter – connects computer to RS-232C port IEEE-488 interface cable connects computer to IEEE-488 interface – cable is 1 m (3.3 ft) long Instrument recalibration with certificate

Lake Shore Cryotronics, Inc.

(614) 891-2244

fax: (614) 818-1600